

Product Specification Sheet

ORDERING INFORMATION

Product name	PetXmark™
Item number	21110590
Product code	2111
IFU languages	English
Sales unit	Box w/ 5 x 1 ml
Box size / total weight	76.5 x 17 x 87.5 mm / 25.5 g



PRODUCT INFORMATION

Indications	PetXmark™ is indicated for use to radiographically mark subcutaneous tissue in dogs and cats.
Product description	<ul style="list-style-type: none"> PetXmark™ is a single use, implantable medical device consisting of a sterile transparent liquid in a one-point-cut (OPC) glass ampoule Each ampoule contains 1 ml liquid (nominal volume) Sterilized using steam Magnetic Resonance (MR) safe Upon injection of the PetXmark™ liquid into soft tissue, efflux of ethanol leads to the formation of a radiopaque, sticky and gel-like fiducial marker <i>in-vivo</i>
Product components	The liquid is a mixture of ethanol, sucrose acetate isobutyrate (SAIB) and an iodinated and acylated derivative of sucrose (x-SAIB).
Storage instructions	Store at 2–25 °C. Keep away from sunlight.
Shelf-life	3 years.
Contraindications	Do not use in patients with known hypersensitivity to iodine or any other component in PetXmark™.
Precautions	<p>PetXmark™ must be handled only by qualified health care professionals. Injection of PetXmark™ requires experience with implantation of fiducial markers or with taking biopsies from the tissue intended for the injection of PetXmark™.</p> <p>See Instruction For Use for detailed precautions and warnings at www.nanovi.com</p>
Appearance	Upon injection of the PetXmark™ liquid into soft tissue, efflux of ethanol leads to the formation of a radiopaque, sticky and gel-like fiducial marker <i>in-vivo</i> .
Usage	<p>PetXmark™ is used as a fiducial marker and is visible on all relevant imaging modalities including MRI.</p> <p>See Instruction For Use for detailed implantation procedure at www.nanovi.com</p>
Compatible equipment	Can be implanted with standard equipment for percutaneous injections.

CONTACT

NANOVI A/S, Diplomvej 378, DK-2800 Kgs. Lyngby, Denmark

Phone: +45 38414059, E-mail: sales@nanovi.com, VAT no. (CVR): 34889031

For more information of this product please visit us at www.nanovi.com